Quick Guide

Data logger for continuous monitoring systems

VDL200

VAISALA

PUBLISHED BY

Vaisala Oyj
Vanha Nurmijärventie 21, FI-01670 Vantaa, Finland
P.O. Box 26, FI-00421 Helsinki, Finland
+358 9 8949 1
vaisala.com
docs.vaisala.com

© Vaisala 2025

No part of this document may be reproduced, published, or publicly displayed in any form or by any means, electronic or mechanical (including photocopying), nor may its contents be modified, translated, adapted, sold, or disclosed to a third party without prior written permission of the copyright holder. Translated document and translated portions of the document are based on the original English version. Parts of the translation may have been created using machine-assisted translation. In ambiguous cases, the English version is applicable, not the translation.

The contents of this document are subject to change without prior notice.

Local rules and regulations applicable to the products and services may vary, and they shall take precedence over the information contained in this document. Vaisala makes no representations on this document's compliance with the local rules and regulations applicable at any given time, and hereby disclaims any and all responsibilities related thereto. You are instructed to confirm the applicability of the local rules and regulations and their effect on the intended use of the products and services.

This document does not create any legally binding obligations for Vaisala towards customers or end users. All legally binding obligations are set forth exclusively in the applicable contract or in the relevant set of General Conditions of Vaisala (vaisala.com/policies).

This product contains software developed by Vaisala or third parties. Use of the software is governed by license terms and conditions included in the applicable contract or, in the absence of separate license terms and conditions, by the General License Conditions of Vaisala Group.

This product may contain open-source software (OSS) components. In the event this product contains OSS components, then such OSS is governed by the terms and conditions of the applicable OSS licenses, and you are bound by the terms and conditions of such licenses in connection with your use and distribution of the OSS in this product. Applicable OSS licenses are included in the product itself or provided to you on any other applicable media, depending on each individual product and the product items delivered to you.

Table of contents

English.....5

Español.....31

Introduction to VDL200 data logger

This document is a quick guide for installation of the VDL200 Data Logger. For the complete user guide, see [VDL200 Data Logger User Guide \(M212938EN\)](#) available at docs.vaisala.com.

VDL200 is an Ethernet-connected data logger for continuous monitoring systems. It provides high-accuracy measurements from up to 2 detachable humidity, temperature, and carbon dioxide probes. VDL200 also supports the ANP115 Analog Input Probe for recording and display of measurements from analog and boolean signals. Suitable applications include real-time measurements in environments where wired connections are preferred. VDL200 is compatible with the Vaisala viewLinc Enterprise Server software and Jade Smart Cloud monitoring solution.


Figure 1 VDL200 and its network connections

VDL200 is powered from the Ethernet connector using Power over Ethernet (PoE). If the network does not provide PoE, a separate PoE power supply must be used. Single-use 1.5 V AA size batteries are used as a temporary backup power source when PoE power is not available. Operating time on battery power depends on the inserted battery type and power usage of the connected measurement probes.

VDL200 maintains accurate time by synchronizing with a Network Time Protocol (NTP) server. VDL200 must be able to synchronize with an NTP server, otherwise it cannot record or send measurement data.

When operating on battery power, the network connection is disabled to save power. If VDL200 has already synchronized its internal clock with a time server, it will keep recording measurement data to its local memory. When PoE power is reconnected, the network connection is restored and the connected monitoring system will automatically retrieve any locally recorded data that is needed.

The following ports are used in the communication between VDL200 and the servers:

- TCP port 8883 (viewLinc Enterprise Server)
- TCP port 443 (Jade Smart Cloud)
- UDP port 123 (NTP servers)

Parts


Figure 2 Front parts

- 1 Display
- 2 Status LED
- 3 Left button
- 4 Right button


Figure 3 Connectors

- 1 Ethernet port (RJ-45)
- 2 Service port (USB-C)
- 3 Probe port 1 (4-pin M8 female)
- 4 Place for locking screw (3×8 mm screw recommended)
- 5 Probe port 2 (4-pin M8 female)


Figure 4 Rear parts

- 1 Power button
- 2 Backup batteries (4 pcs 1.5 V alkaline (LR6) or lithium (FR6))
- 3 Battery cover
- 4 Type label area


Figure 5 Mounting bracket parts

- 1 Screw mounting holes
- 2 Cable tie mounting holes
- 3 Mounting magnet

Display


Figure 6 Display layout

- 1 Number of the probe port that is the source of this measurement
- 2 Abbreviation of measurement parameter
- 3 Latest measured value
- 4 Unit of measurement
- 5 Information icon. Touching the left button (below the icon) opens the device information page.
- 6 Service port connection status
- 7 PoE power supply status
- 8 Host system connection status
- 9 Battery level


Figure 7 Display when errors are active

- 1 Error icon. When the error icon is shown, touching the left button opens the error page.
- 2 PoE power supply status. In this example PoE is not connected.
- 3 Host system connection status. In this example the Ethernet connection is not available as PoE is not connected.
- 4 The color of the icon row is inverted when there are active errors.

Specifications

Table 1 VDL200-compatible probes

Measurement or function	Probe models
Humidity and temperature (RH+T)	HMP110, HMP115
Temperature (T)	HMP110T, HMP115T, TMP115
Carbon dioxide (CO ₂)	GMP251, GMP252
Analog input and door switch	ANP115, ANP115C

Table 2 VDL200 operating environment

Property	Description/Value
Operating environment	Indoor use
Use in wet location	No
Storage temperature	-30 ... +60 °C (-22 ... +140 °F)
Operating temperature	-20 ... +60 °C (-4 ... +140 °F)
Operating humidity	0-100 %RH, non-condensing
IP rating	IP30: Protected against solid foreign objects of 2.5 mm Ø and greater.
Do not place the VDL200 unit in an environment that can exceed this specification, for example inside a climate chamber. Insert only the measurement probe(s) in the chamber and leave the VDL200 unit outside it.	

Table 3 VDL200 powering

Property	Description/Value
Powering options	Power over Ethernet (PoE) Battery backup ¹⁾ USB-C for temporary service use
PoE supply voltage	48 V DC
USB-C supply voltage	5 V DC
Power consumption	Max. 2 W
Battery type	AA size, 1.5 V, LR6 (alkaline) or FR6 (lithium)
Number of batteries	4
Typical operating time on battery power at +20 °C (+68 °F)	
RH and T measurement combinations	2 weeks with alkaline batteries
CO ₂ measurement or analog input	24 hours with lithium batteries

1) Network connection is not available when VDL200 is operating on battery power.

Table 4 VDL200 inputs and outputs

Property	Description/value
Ethernet connector	8P8C (RJ-45)
Probe connector (2 pcs)	4-pin female M8 connector
Service port	USB-C
Ethernet interface	
Supported standards	10BASE-T, 100BASE-TX
IPv4 address assignment	DHCP (automatic), static IP

Status icons

Table 5 Status icons on VDL200 display

Icon(s)	Icon name	Meaning
	Information	Access the device information page by touching the button below this icon.
	Error	A device error is active. Touch the button below this icon to access the error page and the device information page.
	Service port active	Service port is in use.
	PoE available	Device is currently powered by Power over Ethernet (PoE).
	PoE not available	PoE is not available and the device is operating on backup battery power. This state is considered an error state. Ethernet connection is disabled to reduce power consumption, which prevents communication with the monitoring system and the time servers.
	Connection OK	Monitoring system address is configured and the system is reachable. This icon is used for viewLinc Enterprise Server.
	Connection OK	Monitoring system is reachable. This icon is used for cloud-based monitoring systems such as Jade Smart Cloud.
	Connection unavailable	Ethernet connection is not available.
	Check connection	Connection to monitoring system not available.
	Battery charge	Approximate charge remaining in the backup batteries.
	Batteries empty	Backup batteries do not have sufficient charge to power the device.
	Batteries missing	Backup batteries are missing or inserted incorrectly.

Error messages

When VDL200 has one or more error messages to show, it flashes the status LED with red, and displays the  icon over the left button. Touching the button opens the **Errors** screen where you can view the error messages.

Some errors are associated with status icons that are shown on the home screen.

Table 6 VDL200 error messages

Error message and icon(s)	Description	Recommended action
No probes detected. Connect probes and restart.	No compatible measurement probes were detected when VDL200 was started up. This message is shown on the home screen instead of the Errors screen.	Connect measurement probes and restart the VDL200.
Power over Ethernet not available. 	PoE power is not available. VDL200 is operating on battery power or power from the temporary service port connection. This error means VDL200 cannot operate its network connection, which causes additional errors to be active.	1. Verify that the Ethernet cable is connected. 2. Enable PoE from the network switch or connect a PoE power source between the network and the VDL200.
Network connection not available. 	VDL200 does not have a network connection.	1. Verify that the Ethernet cable is connected. 2. Connect the VDL200 to Insight software and verify the network settings. 3. Depending on how the network is configured, you may need to supply your IT support with the MAC address of the VDL200 so that it can be allowed to the network.
Monitoring system connection problem. 	If VDL200 is used with a cloud-based monitoring system, this error is likely caused by lack of Internet connectivity.	There may be a connectivity problem between the monitoring system and VDL200. Request help from your local IT support to resolve it.
	If VDL200 is used with viewLinc Enterprise Server, the address of the monitoring system must be manually configured. The address may be missing or incorrect.	Connect VDL200 to Insight software and verify that the address of the monitoring system has been configured correctly.

Error message and icon(s)	Description	Recommended action
Backup batteries empty. 	The battery level is insufficient for operation.	Replace the batteries.
Backup battery level low. 	The battery level is low but currently sufficient for operation.	
Backup batteries not installed. 	Backup batteries are not installed.	Insert new batteries.
	Backup batteries are installed in incorrect orientation.	Verify the correct orientation of the batteries from the markings on the battery compartment, and correct the orientation.
	The battery isolation tab that prevents battery drain during storage has not been removed.	Remove the battery isolation tab.
System time not available.	The real-time clock of the VDL200 does not have system time set. This error means the VDL200 cannot store and send measurement data.	Resolve the time synchronization error.
Time synchronization problem.	VDL200 is currently unable to synchronize its system time with any of the network time protocol (NTP) servers.	1. Verify that VDL200 has a network connection and PoE power connected. 2. Connect the VDL200 to Insight software and verify that the addresses of the NTP servers have been configured correctly. 3. A network firewall may be blocking UDP port 123 that is needed for the NTP protocol to function.

Installing VDL200

To install VDL200, you need the following tools and equipment:

- Content of the VDL200 product package:
 - VDL200 device with mounting bracket
 - According to order: measurement probes, probe cables, probe holders, and batteries
- If not included in the order: 4 pcs AA size 1.5 V batteries, type LR6 (alkaline) or FR6 (lithium)
- Cable ties or screws and wall plugs for mounting
- Cable tie cutter or screwdriver, depending on mounting method
- Ethernet cable
- Power over Ethernet (PoE) power source

For configuring the settings of VDL200:

- Computer with Microsoft Windows® operating system (64-bit version) and Insight PC software (version 1.3.0 or later) installed, and a free USB port
- USB-C cable for connecting VDL200 to a computer


NOTICE! VDL200 must synchronize its real-time clock with a time server before it can store and send data. As long as the error **System time not available** is active, no new measurement data is stored or sent to the monitoring system. Complete the installation steps to remove the error.

- ▶ 1. If you are using the ANP115 Analog Input Probe to convert an analog measurement signal, configure it using Insight PC Software before connecting it. See ANP115 documentation at vaisala.com/anp115/docs.

The door switch model ANP115C is preconfigured and ready to use.

2. Slide the mounting bracket out of the data logger.


WARNING! Magnetic hazard

The mounting bracket has a strong magnet. Magnets can cause serious harm by interfering with the operation of pacemakers and other implanted medical devices. Keep the magnet away from persons with such devices, and from any other devices that are sensitive to magnetic fields. Do not separate the magnet from the mounting bracket.

3. Connect the measurement probes to the probe ports of the data logger:

- Direct connection (ANP115/C, HMP115/T, and TMP115 only): Verify the orientation of the connector pins and push the probe to the probe port. Secure the probe using a probe nut.
- Cable connection (any probe): Connect the cable to the probe port, and secure it using the locking ring of the cable connector. Then connect the probe to the cable.

4. If you are using the ANP115 Analog Input Probe, connect its input cable to the source of the analog measurement signal. If using the ANP115C, mount the door switch in the intended location. Refer to ANP115 documentation for wiring and connection instructions.
5. Open the battery cover of VDL200.

6. If the batteries have been inserted at the factory, verify their type (LR6 alkaline or FR6 lithium) and remove the battery isolation tab if it is present. The tab prevents the batteries from draining during storage.
7. If the battery compartment is empty, insert 4 pcs AA size 1.5 V batteries, type LR6 (alkaline) or FR6 (lithium). Follow the markings on the battery compartment to insert the batteries in the right orientation.

**WARNING! Explosion and combustion hazard**

Use of incompatible batteries, inserting batteries incorrectly, and exposing batteries to high temperatures may cause batteries to explode, deflagrate, or leak flammable liquid or gas, which may lead to serious injury. Always verify the required battery type and the correct orientation inside the battery compartment. Also, make sure the temperature at the installation location is within the operating temperature range of the batteries.

Use batteries that have a full charge and are not past their expiration date. Lithium batteries must be used with CO₂ measurement, and when the ambient temperature of the data logger's installation location is below 0 °C (+32 °F).

8. Press the power button for approximately 1 s to turn on the data logger. The data logger plays a short sound effect when turning on and off.

The data logger starts up on battery power and detects the connected probes. The measurements from the probes are automatically assigned for display and recording.


Batteries are intended for backup use only. Proceed to connect Power over Ethernet (PoE) as soon as possible. **The Ethernet connection is disabled until PoE is available.**

9. Replace the battery cover.
10. Connect the Ethernet cable.
11. If the network does not provide PoE, connect a PoE power source between the network and the data logger.
12. Connect the VDL200 to your laptop using a USB-C cable and configure its settings. See [Configuring settings using Insight software \(page 18\)](#).

13. Attach the mounting bracket to the installation location with screws, cable ties, or the magnet.

Note the orientation – when the data logger is installed, the probe ports should point down.

14. Slide the data logger into the mounting bracket.
15. Optional: Insert a 3×8 mm screw into the locking screw hole to prevent the data logger from being easily removed.
16. Use the included accessories to attach cabled probes to their final locations and secure the cables. See [Mounting probes \(page 20\)](#).
17. Log in to your monitoring system to as an administrator to verify that the VDL200 is connected and finalize device configuration. See the applicable section for your monitoring system:
 - [Remote management with viewLinc Enterprise Server \(page 26\)](#)
 - [Remote management with Jade Smart Cloud \(page 27\)](#)

Remote management provides additional configuration options that may be relevant for you. For example, the temperature unit displayed by the data logger can be changed remotely. The default temperature unit is degrees Celsius.

18. Check the display of the VDL200 to verify there are no active errors, and that the expected measurements are displayed.

If an ANP115 Analog Input Probe is connected, verify that the device that sends the analog output signal is powered on, and the selected parameter and unit are shown with a valid value. Note that the display uses abbreviations for parameter names.

For the list of abbreviations, see [VDL200 Data Logger User Guide \(M212938EN\)](#).

If the door switch model ANP115C is connected, operate the door or lid that is monitored by the switch to verify that the value on the display changes.

Configuring settings using Insight software

Required tools:

- Computer with Microsoft Windows® operating system (64-bit version) and Insight PC software (version 1.3.0 or later) installed, and a free USB port
- USB-C cable for connecting VDL200 to a computer


Figure 8 Configuring VDL200 settings with Insight

1. Turn on the VDL200 if it is off.
2. Open the Insight software on your computer.
3. Connect the USB cable between your computer and the service port of the VDL200.
4. Wait for the Insight software to detect the data logger.
5. In Insight, select  > **Configure device**.
6. Select the **Battery type** so that it matches the batteries currently inserted in the device.

The battery type setting is important to optimize the accuracy of the battery indicator and to allow the monitoring system to reliably send the low battery warning.

7. Configure the **Network** setting to select addressing type:

- To use automatic IP address configuration using DHCP, select **Automatic (DHCP)**.
- To configure a static IP address, select **Static IP**.

Your network administrator may request the **MAC address** of the device to allow it to the network. It is marked on the device itself, and you can also press the left button on the device to view it on the device information screen.

8. If you selected Static IP configuration, fill the following fields as instructed by your network administrator:

- **IP address**
- **Subnet mask**
- **Gateway**
- **DNS server 1**

9. Only when used with viewLinc Enterprise Server: Configure the hostname, fully qualified domain name, or IP address of your viewLinc ES system in the **Monitoring system address** field.

VDL200 communicates with viewLinc ES on TCP port 8883 by default. You do not have to specify the default port in the address. However, if the port has been changed on the viewLinc side, you must add a colon ":" and the correct TCP port number to the address. For example:
viewlinc.examplecompany.com:1234

10. If you have an NTP server (time server) in your local network, overwrite **NTP address 1** with its hostname, fully qualified domain name, or IP address.

11. Select **Save** when done.

Mounting probes

For more information on probe mounting and available accessories, refer to the user guides of the probes:

- [HMP60 and HMP110 Series User Guide \(M211060EN\)](#)
- [GMP251 and GMP80P User Guide \(M211799EN\)](#)
- [GMP252 User Guide \(M211897EN\)](#)
- [Vaisala ANP115 Analog Input Probe Quick Guide \(M213202EN\)](#)

Probe holder



Figure 9 Probe holder ASM213382

- 1 Holes for attaching the holder with screws (screws not included)
- 2 Magnet
- 3 Probe attachment

The probe holder is a versatile mounting accessory for securing Ø 12 mm diameter cable mounted probes.

To attach the probe, simply press the body of the probe into the probe attachment.

Figure 10 Example probe in the probe holder

CO₂ probe mounting kit

When the data logger is purchased with a CO₂ measurement probe, a CO₂ probe mounting kit (Vaisala item ASM214253SP) is included. The kit includes a versatile probe support accessory that has attachments for a Ø 25 mm probe (for example, GMP251), a Ø 12 mm probe (for example, HMP110), and the sensor tip of the TMP115 probe. It also provides sufficient separation between the probes to prevent the mild heating effect of the CO₂ probe from affecting the measurement of the other probe. For this reason, it is best not to attach the CO₂ probe to another probe using a cable tie, for example.


Figure 11 Probe support accessory mounting options

- A** Mounting through a hole using the screw-on attachment part. This method can be used to install the probes under a grid shelf, for example.
- B** Mounting using a reusable fastener strip. Clean the attachment surfaces using the included cleaning pad before applying the strip.
- C** Mounting with screws and wall plugs.

Mounting HMP110 probes

HMP110 is a robust stainless steel probe for humidity and temperature measurement in demanding conditions. Suitable for measurement inside chambers, fridges, and freezers in temperature range $-40 \dots +80 \text{ }^{\circ}\text{C}$ ($-40 \dots +176 \text{ }^{\circ}\text{F}$). Must be connected to the data logger using a cable. Probe diameter 12 mm (0.47 in).


Figure 12 HMP110 probe

- 1 M12x1 thread for through-wall installation using mounting nuts.
- 2 Attach from this area using probe holder ASM213382 or cable tie.
- 3 Sensor protection filter. Do not attach from this area.

Mounting HMP115 probes


Figure 13 HMP115 probe

- 1 If using a probe holder, align it to this groove.
- 2 Plastic sleeve that locks the probe in place when connected to a mechanically compatible host device (for example, RFL100 data logger). Diameter 14 mm (0.55 in) at this point.
- 3 Attach from this area if using a cable tie. Diameter 12 mm (0.47 in) at this point.
- 4 Sensor protection filter. Do not attach from this area.

Mounting TMP115 probes


Figure 14 TMP115 probe

- 1 If using a probe holder, align it to this groove.
- 2 Plastic sleeve that locks the probe body in place when connected to a mechanically compatible host device (for example, RFL100 data logger). Diameter 14 mm (0.55 in) at this point.
- 3 Probe body with measurement electronics. Attach from this area if using a cable tie. Diameter 12 mm (0.47 in) at this point. Operating temperature range $-40 \dots +60 \text{ }^{\circ}\text{C}$ ($-40 \dots +140 \text{ }^{\circ}\text{F}$).
- 4 Sensor cable. Do not cut or bend into a tight loop.
- 5 Sensor tip, diameter 4.8 mm (0.19 in). Secure using a cable tie or insert into thermal dampener block for added thermal mass. Insert into the probe support accessory when using together with a CO_2 probe. Operating temperature range $-196 \dots +150 \text{ }^{\circ}\text{C}$ ($-320.8 \dots +302 \text{ }^{\circ}\text{F}$).

CAUTION! The operating temperature range of the sensor tip is much wider than that of the probe body. Leave the probe body outside the measured environment if possible, and avoid inserting it in environments that are outside its operating range.

When working with equipment in extremely cold temperatures, use appropriate personal protective equipment such as thermally insulated gloves and clothing. Wear protective eyewear if working with coolants such as liquid nitrogen, and observe safe handling and storage precautions.

Mounting GMP251 and GMP252 probes

GMP251 and GMP252 are robust carbon dioxide (CO₂) measurement probes for use in demanding applications such as life science incubators. Operating temperature range -40 ... +60 °C (-40 ... +140 °F). The CO₂ probe mounting kit (Vaisala item ASM214253SP) is recommended for mounting, but the probes can also be attached from the probe body by other methods.


Figure 15 GMP251 probe


Figure 16 GMP252 probe

- 1 M12 5-pin male connector. Connect to VDL200 using the included M8-M12 adapter cable (item code CBL211291SP).
- 2 Ø 25 mm probe body.
- 3 Sensor protection filter. Do not attach from this area.

Mounting ANP115 Analog Input Probe

ANP115 converts a current, voltage, or contact input from an analog sensor into a digital output compatible with the VDL200 Data Logger. It supports a wide range of passive and active contacts. You must configure the ANP115 using Vaisala Insight software before installation.

Operating temperature range of ANP115 is -40 ... +60 °C (-40 ... +140 °F).


Figure 17 ANP115 open wires model

- 1 Probe nut for mounting directly to VDL200
- 2 Attach from this area if using a probe holder or cable tie. Diameter 12 mm (0.47 in) at this point.
- 3 Instruction label with wiring options
- 4 Open wires

ANP115C is a door switch model that is preconfigured at the factory. A mounting kit for the door switch is included with the ANP115C. It contains some mounting tape and the magnet that operates the switch.

Operating temperature range of ANP115C is $-20 \dots +60 \text{ }^{\circ}\text{C}$ ($-4 \dots +140 \text{ }^{\circ}\text{F}$).


Figure 18 ANP115C door switch model

- 1 Probe nut for mounting directly to VDL200
- 2 Attach from this area if using a probe holder or cable tie. Diameter 12 mm (0.47 in) at this point.
- 3 Magnetic reed switch

Remote management with viewLinc Enterprise Server

- ▶ 1. Open the viewLinc Enterprise Server user interface by double-clicking the desktop icon, or by entering the server's IP address and port in a web browser's address field.
- 2. Log in with a user account that has the right to manage devices.
- 3. When viewLinc detects new devices, **New Devices** text appears at the top of the screen. If there are any new devices:
 - a. Click the **New Devices** text to open the **New Devices** window.
 - b. Compare the list of new devices to your device installation plan to verify that all your devices have been found. In the **Accept** column, select each device you want to add.
 - c. Select **Save** to register the selected devices with your viewLinc Enterprise Server.
- 4. To access management options for accepted devices:
 - a. Select **Sites Manager > Hosts and Devices**.
 - b. Locate the device you want to manage from the device tree, right-click the entry of the device, and select **Edit Properties**. The **Edit Device Properties** window opens.
 - c. In this window you can set some information fields that viewLinc uses (for example, **Device alias**) and some actual settings that are sent to the device. Note the following settings for VDL200:
 - **Metric units enabled**: Change the setting to **No** to display degrees Fahrenheit for temperature on the device display.
 - **Device log level**: Do not adjust this setting unless instructed by Vaisala to do so. Changing the setting will increase the power consumption of the device, which will reduce the battery backup time.

Remote management with Jade Smart Cloud

- ▶ 1. Log in to your Jade Smart Cloud account at jade.vaisala.com.
- 2. Select the **Devices** menu and locate the device you want to manage on the list.
- 3. Select the device from the list. The **Overview** tab opens, showing the communication status of the device.
- 4. Select the **Device properties** tab.
- 5. Check the **Units on display** value. This value controls the unit system that is used on the local display. To change the setting:
 - a. Select **Edit**.
 - b. Select the new setting from the drop-down menu.
 - c. Select the check mark to submit the change.

Vaisala services

Service offering and support channels

See our service offering at vaisala.com/services.

See our support channels at vaisala.com/support.

Service centers

Vaisala offers comprehensive customer care throughout the product life cycle. Our service centers offer maintenance, calibration, upgrade, and repair services. All the services are provided worldwide with fast deliveries. For more information, go to vaisala.com/services.

The services are available from:

- Vaisala Online Store, store.vaisala.com
- Vaisala Services, helpdesk@vaisala.com

Technical support

You can contact Vaisala technical support at helpdesk@vaisala.com. Provide the following supporting information as applicable:

- Product name, model, and serial number
- Software/Firmware version
- Name and location of the installation site
- Name and contact information of a technical person who can provide further information

When contacting technical support for troubleshooting, also provide the following information as applicable:

- What failed (what worked / did not work)?
- Where did it fail (location and environment)?
- When did it fail (date, immediately / after a while / periodically / randomly)?
- How many failed (only one defect / other same or similar defects / several failures in one unit)?
- What was done when the failure was noticed?
- What was connected to the product and to which connectors?
- What is the input power source type, voltage, and what else (such as lighting, heaters, and motors) was connected to the same power output?
- Are all parts connected and grounded properly? Take a photo to help the troubleshooting.

If available, provide the diagnostics package for the application.

Warranty

For warranty terms and conditions and the warranty claim process, see vaisala.com/warranty.

Note that warranty coverage may not apply if the product is damaged due to normal wear and tear, exceptional operating conditions, negligent handling or installation, or unauthorized modifications. See the applicable supply contract or Conditions of Sale for details of the warranty for each product.

Manufacturer and importer information

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Name of equipment: VDL200 Data Logger

Country of manufacture: Finland

Link to RRA registration information: www.rra.go.kr/selform/rws-VDL200

FCC Part 15 compliance statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that the interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

CAUTION! This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

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Kenelec Scientific Pty Ltd