



## Industrial Tomography Systems extends the University discount Program to the U2s Ultrasound spectroscopy system

Since 2009, Industrial Tomography Systems (ITS) through their University Discount Program have offered University Researchers affordable and advanced in-situ Analytical Instrumentation. This program has now been expanded to include ITS' newly launched U2s Ultrasound Spectroscopy System. This cutting edge technology will allow researchers to:

- **Convert Ultrasound Spectra to useful Viscosity Data in real time**
- **Determine viscosity of samples in-situ – either in a vessel or in a pipe**
- **Characterize rheology of emulsions, slurries and soft solids**
- **Perform in-situ measurements with no limitation to the sample opacity**
- **Measure changes to viscosity rates of up to 10 spectra/second**



### University Discount Offering:

The U2s with standard software, flow cell sensor and/or probe, full installation and two days of training will be offered to Universities at a significant discount.

Further discounts are available for multiple purchases.

**To schedule a free demonstration to characterize your samples at your University:**

Please call or e-mail:	<b>Paul Montague</b> (02) 8977 4017 paul@kenelec.com.au
------------------------	---