



The TP-Series – the economical line for the Lab



The new TP-series – simply easy, simply good, simply economic

A new range of compact balances designed to lead the way in providing precision weighing at a great price. Extensive application programs, well-established weighing technology and easy handling makes these balances an ideal solution for standard weighing tasks in the lab, field, industrial or educational environments. Quick accurate weighing powered by state-of-the-art microprocessor technology ensures balance stability and reproducible results.

- Large easy to read high contrast LCD
- Sealed easy to clean keypad
- Stainless Steel weighing pan
- Compact and sturdy design
- Level Indicator and adjustable feet
- Internal, manual calibration in the TP-303i
- Automatic calibration with external weight
- Applications include unit conversion, parts counting, % weighing, display hold, totalizing
- Capable of displaying 14 multiple modes
- Security eyelet
- Bi-directional RS-232 interface
- Made in Germany
- 3 years warranty
- Accessories: Printer, Cal-weights, Dust cover

| Model | Capacity | Readability | Pan size | Stability time |
|---------|----------|-------------|--------------|----------------|
| TP-214 | 210 g | 0.1 mg | Ø 80 mm | 3 sec |
| TP-303i | 310 g | 0.001 g | Ø 97 mm | 3 sec |
| TP-303 | 310 g | 0.001 g | Ø 97 mm | 3 sec |
| TP-3002 | 3100 g | 0.01 g | 174 × 143 mm | 2.5 sec |
| TP-1502 | 1500 g | 0.01 g | 174 × 143 mm | 2.5 sec |
| TP-6101 | 6100 g | 0.1 g | 174 × 143 mm | 2 sec |
| TP-3101 | 3100 g | 0.1 g | 174 × 143 mm | 2 sec |
| TP-12 | 12000 g | 1 g | 174 × 143 mm | 1.5 sec |

Following balances are supplied as EU-stamp approval ex factory

| | | | | |
|----------|--------|--------|--------------|-------|
| TP-214A | 210 g | 0.1 mg | Ø 80 mm | 3 s |
| TP-124A | 120 g | 0.1 mg | Ø 80 mm | 3 s |
| TP-512A | 510 g | 0.01 g | Ø 116 mm | 2 s |
| TP-1200A | 1200 g | 1 g | 174 × 143 mm | 1.5 s |



TP-214, TP-214A, TP-124A



TP-303



TP-303i



TP-3002, TP-1502, TP-6101,
TP-3101, TP-12

UK office:

Denver House, Sovereign Way
Trafalgar Business Park
Downham Market, Norfolk PE38 9SW
Phone +44.1366.386-242
Fax +44.1366.386-204
uk.sales@denverinstrument.com

Denver Instrument GmbH

Robert-Bosch-Breite 10
37079 Goettingen, Germany
Phone +49.551.209 77-30
Fax +49.551.209 77-39
www.denverinstrument.com
info@denverinstrument.de