

The Low-Budget Series PK and XP Compact Balances of the Highest Quality



The Low-Budget Series – Precision Balances

These balances are always within reach when you need them. In the lab and in the field, for quality control work and use at universities, schools and jewellery stores – the low-budget series from Denver Instrument is the right choice for all environments in which precision weighing is required at an affordable price.

- Quick, accurate weighing
- State-of-the-art microprocessor technology ensures
- Balance stability and reproducible results
- Keys “click” audibly when pressed, making for easy and sure operation

High-quality features

- High-contrast, large, easily readable LCD display
- Sealed, easy-to-clean keypad with large tactile-response keys
- Stainless steel weighing surface
- Rugged housing made of sturdy ABS plastic
- Compact design
- Two-year warranty

Model	Readability	Weighing range	Pan size
PK-352	0.01 g	350 g	ø 115 mm
PK-202	0.01 g	200 g	ø 115 mm
PK-4801	0.1 g	4,800 g	141 x 145 mm
PK-2401	0.1 g	2,400 g	141 x 145 mm
PK-1201	0.1 g	1,200 g	141 x 145 mm
PK-601	0.1 g	600 g	141 x 145 mm
PK-401	0.1 g	400 g	141 x 145 mm
PK-10	1 g	10 kg	141 x 145 mm
PK-6	1 g	6 kg	141 x 145 mm
PK-4	1 g	4 kg	141 x 145 mm
XP-300	0.01 g	300 g	127 x 178 mm
XP-600	0.02 g	600 g	127 x 178 mm
XP-1500	0.05 g	1,500 g	127 x 178 mm
XP-3000	0.1 g	3,000 g	127 x 178 mm



PK 4 – 4801

PK 352/202

XP 300 – 3000

for all Weighing Needs

PK Series

- 9V block battery
- 4 weighing modes: gram, ounces, carat, pennyweight
- Counting program, formulation program and memory function
- Automatic calibration with external weight (weight for PK-352/202/401/601 included)
- Level indicator and adjustable feet
- Security eyelet for theft-proof mounting
- Adapter as a standard feature, 9V block battery optional (except PK 352/4801/2401)
- Optional RS232C interface
- Overload protection

XP Series

- Integrated rechargeable battery
- 7 weighing modes: gram, ounces, troy ounces, pennyweight, carat, grain or tael
- Automatic calibration
- Part counting
- Percentage weighing
- Overload protection
- Battery-saving sleep mode
- Optional RS232C interface



Denver Instrument GmbH
Robert-Bosch-Breite 10
37079 Göttingen
Phone +49.551.20977-30
Fax +49.551.20977-39
www.denverinstrument.com
info@denverinstrument.de

UK office:
Denver House, Sovereign Way
Trafalgar Business Park
Downham Market, Norfolk PE38 9SW
Phone +44.1366.386-242
Fax +44.1366.386-204

Not only for laboratory use:
with its counting program, the
balance also cuts a fine figure
in the stockroom.

