

# DLT – Series

Best value for money



## DLT – Series Best value for money

Denver Instrument offers with the DLT-line a new series of portable, compact and very economically priced balances for all basic weighing tasks in this segment. The DLT-line is a real new alternative for all demands in schools, in the Lab or in the factory, where simple and cheap weighing devices are required.

Despite the very attractive prices the Denver DLT offers very reliable and precise weighing results and easy handling with just 3 buttons in the well-known Denver tradition.

### Features and Equipment

- Large display
- Sturdy housing and compact design
- Battery – or AC powered (power supply optional)
- Energy saving mode
- Easy handling
- External Calibration
- Balances are stackable
- 9V battery included
- Antitheft locking device
- Battery saving mode
- Transportation locking device
- 2 Years warranty



Model	Capacity	Read-ability	Measuring time	Linearity	Pan size	Netweight
DLT-411	410 g	0.1 g	1.5 sec	± 0.2 g	135 × 140 mm	440 g
DLT-211	210 g	0.1 g	1.5 sec	± 0.2 g	135 × 140 mm	440 g
DLT-4100	4,100 g	1 g	1.5 sec	± 2 g	135 × 140 mm	440 g
DLT-2100	2,100 g	1 g	1.5 sec	± 2 g	135 × 140 mm	440 g

Denver Instrument GmbH  
Robert-Bosch-Breite 10  
37079 Goettingen, Germany  
Phone +49.551.209 77-30  
Fax +49.551.209 77-39  
www.denverinstrument.com  
info@denverinstrument.de

Denver House, Sovereign Way  
Trafalgar Business Park  
Downham Market, Norfolk PE38 9SW  
Phone +44.1366.386-242  
Fax +44.1366.386-204  
www.denverinstrument.com  
info@denverinstrument.de