

ENGINE EXHAUST CONDENSATION PARTICLE COUNTER MODEL 3790

STANDARD REPLACEMENT PARTS

Standard replacement parts for Model 3790 Condensation Particle Counter consist of:

<u>Part Number</u>	<u>Part Description</u>
390065	Aerosol Instrument Manager® Software for CPC/EAD
1031515	Maintenance Kit for 3772 and 3771 (includes 1031513 and 1031514)
1031513	Replacement Filter Kit for 3772 and 3771 (includes 2 micropump filters and 3 butanol fill/drain filters)
1031514	Replacement Saturator Wick Kit for 3772 and 3771 (qty. 2 saturator wicks)
1500108	SD/MMC Data Memory Card
1030008	3771/3772 Critical Orifice Assembly
1030010	3771/3772/3010/3760A Replacement Nozzle Assembly
1031548	3771/3772 CPC Fill/Drain Replacement Kit CPC (includes fill bottle, drain bottle, and bottle bracket)
1031492	Large Charcoal Filter Kit for CPCs (includes 5 large charcoal filters)
1031493	Small Charcoal Filter Kit for CPCs (includes 5 small charcoal filters)

For pricing information, or to place an order you can 1) send the necessary information to answers@tsi.com or 2) contact your local TSI representative (see contact information below) and speak to one of our customer specialists. Please include your contact information (telephone number and mailing address) with all requests.

Aerosol Instrument Manager® is a registered trademark of TSI Incorporated.



UNDERSTANDING, ACCELERATED

TSI Incorporated – Visit our website www.tsi.com for more information.

USA	Tel: +1 800 874 2811	India	Tel: +91 80 67877200
UK	Tel: +44 149 4 459200	China	Tel: +86 10 8251 6588
France	Tel: +33 4 91 11 87 64	Singapore	Tel: +65 6595 6388
Germany	Tel: +49 241 523030		