

# MITA MASK ADAPTER KIT MODEL 8125-06

---

INSTRUCTION GUIDE



APPLICABLE MASK MODEL:  
3M 6000 SERIES, 7000 SERIES, 7800S, 7800B

## Included Accessories

**40/74 mm Mask Filter Plug – Qty 2**  
(for 3M 7000 series w/40mm threaded connections)

- Installed on mask during **Protection Factor** test



**Mask Filter Plug – Qty 2**  
(for 3M masks w/bayonet connections)

- Installed on mask during **Protection Factor** test



# Components Installed on Masks



Mask Filter Plug

Example of 3M 6900



Mask Filter Plug

Example of 3M 7800



40/74 mm  
Mask Filter Plug

**Example of 3M 7800**



UNDERSTANDING, ACCELERATED

**TSI Incorporated** – Visit our website [www.tsi.com](http://www.tsi.com) for more information.

**USA** Tel: +1 800 874 2811  
**UK** Tel: +44 149 4 459200  
**France** Tel: +33 4 91 11 87 64  
**Germany** Tel: +49 241 523030

**India** Tel: +91 80 67877200  
**China** Tel: +86 10 8251 6588  
**Singapore** Tel: +65 6595 6388