


Follow the below steps to download software for your specific unit from the TSI website.

## Step 1

Go to the TSI website: <https://tsi.com/>

## Step 2

In the header menu, hover over **Software**, then in the drop-down menu click **TSI Software and Firmware Wizard**:




## Step 3

In the text box, **begin typing unit details** then **select the applicable make / model** from the drop-down menu:

PortaCount™ Respirator Fit Tester 8048 [8048]  
PortaCount™ Respirator Fit Tester 8040 [8040]

## Step 4

Click the **Search** button:



## Step 5

Locate the relevant software and **click on the button** to initiate the download process:

### FitPro™ Ultra 6.0.1 Software Setup & FPUClient 4.12.0 Setup



Complete install including .Net support, etc. For Windows® 10 (32-bit & 64-bit) or Windows 11 operating system. FitPro Ultra is not supported on Windows for the ARM64 architecture.

Follow any further instructions provided by the TSI website.